

Table SI. Conditions used for cell irradiation in 25 cm² culture flasks containing 3 ml of cell culture medium.

Activity (Bq)	Irradiation dose (Gy)	Area (cm ²)
1,1x10 ⁴	0,35	25
5,2x10 ⁴	1,7	
1,1x10 ³	3,5	
2,0x10 ⁵	6,6	
3,2x10 ⁵	10	
6,3x10 ⁵	20	
9,3x10 ⁵	30	
1,9x10 ⁶	60	

Table SII. Conditions used for cell irradiation in 75 cm² culture flasks containing 7 ml of cell culture medium.

Activity (Bq)	Irradiation dose (Gy)	Area (cm ²)
7,4x10 ⁴	1	75
1,9x10 ⁵	5	
7,4x10 ⁵	10	
1,5x10 ⁶	20	
4,3x10 ⁶	60	